User manual
Benutzerhandbuch
Manual de usuario
Gebruikershandleiding

Microscopes

Article code | Artikelnummer | Código del artículo | Artikelnummer:

2760611

BMS 100 FL LED


WWW.BMSMICROSCOPES.COM

BMS_100FL_serie_GB-DE-ES-NL_manual_v1. PDF9001100

GB

Thank you for choosing BMS Microscopes.

The BMS 100 FL microscope makes it possible to study structures and textures of botanical and zoological nature at cell level. Both prepared slides or self-made slides can be observed. The material used should allow light to pass easily, too thick material will result in a blurred and undefined image. BMS 100 FL is developed according to the latest technology on optical, electrical and mechanical details.

Including: dust cover, manual and tension wrench

Package dimensions (Width x Height x Depth): 26 x 46 x 18 cm

Gross weight: 2,9 kg

Warranty: 5 years

GETTING STARTED

The objectives have already been placed in the nose piece in the right sequence. They are parfocal and parcentered. Unless necessary, do not remove the objectives. Insert the plug into an earthed power point (check if the voltage matches the microscope). Switch on the light. The microscope is now ready for use.

POWER OPTIONS

- Main power cable (100...240V~/50Hz), available for EU (item no. 76995), UK (item no. 76999), CH (item no. 76998) and USA (item no. 76996).
- Rechargeable batteries (AA, 1,5V), optional item, can be ordered as a pack of three (item no. 29554). When this battery pack is placed, it will be charged automatically when the microscope is connected to the mains and the power switch of the microscope is switched on. When the LED indicator, located on the back of the microscope, is green, the battery pack is fully charged. The microscope can now be used, disconnected from the mains, for 10 hours. When the LED indicator is red, the batteries need to be charged again. Charging time is about 3 hours.
- Non rechargeable batteries (AA 1,5V, 3 pieces). The microscope will also work with these batteries. **WARNING:** do not connect the microscope to the main electricity when these batteries are placed. This will cause severe damage to the microscope.

DOUBLE SAFETY IS GUARANTEED

- Full separation from the main power through DPDT switch.
- Indicator LED to show if main power is switched on and if optional battery is charged.
- Potentiometer for regulation of the High Power LED.
- High Power LED for bright image.
- DPDT (double pole, double throw) when switched off, the power is completely shut off. The BMS 100 FL microscope has a LED illumination light source. This light source is characterized by low maintenance and low power consumption. This LED source has an average life span of 50,000 hours. Besides that, this LED illumination system lacks heat development and is therefore safe to be used by children. It is recommended though that an adult is supervising at all time.

SAFETY INFORMATION

- Think of the safety of yourself and those being around you.
- Always disconnect mains plug from power point before maintenance.
- Work in a clean and proper environment.
- Clean all parts coming in contact with the observed specimen to prevent infection.
- Place the microscope on a sturdy, level table or working surface.
- Always use original accessories.

CAUTION

The delicate components may get damaged if the microscope is not handled with care. Use the microscope in a clean environment. Do not use it when it is subject to direct sunlight, high temperatures, humidity or dust. Carry the microscope with care, do not hold it by the stage. Remove specimen before transport.

REPLACING FUSE

Always disconnect mains plug from the power point before any kind of maintenance. The BMS 100 FL microscope is equipped with one glass fuse, 1,0 A (applicable for item 2760611)

- Wait until the microscope and electrical parts have cooled down.
- The fuse is located at the back of the microscope.
- Open the housing of the fuse by turning it counter clockwise.
- Take broken fuse out.
- Replace new fuse after being sure of using same rating in new fuse.
- Screw fuse holder clockwise back in place.

CARE & MAINTENANCE

- Microscope is a precision instrument, it should be handled with care.
- The microscope should not be in direct sunlight. But placed in a dry and clean place.
- Avoid high temperature and vibration.
- Avoid dirt and fingerprints on the lens surface, this way your image will not be affected.
- All lenses must be kept clean. Fine dust on the surface of the lens should be blown off with hand blower or wiped off gently with a soft lens tissue. Fingerprints or oil mark on it should be wiped off with a tissue moistened with a small amount of xylene or a mixture of alcohol and ether.
- Never use the organic solution to clean the other surface (especially the plastic surfaces). If necessary, please choose the neutral detergent.
- All lenses are mechanically structured and adjusted to the best state for the microscope, please do not disassemble.
- The rack and pinion should be regularly greased with a small amount.
- After use, cover the microscope with the provided dust cover and store it in a dry and clean place, that is free from moisture to prevent rust.

SPECIFICATIONS


For more specifications, please view the product page on the website: www.bmsmicroscopes.com

It is within our company policy to develop continuously. BMS Microscopes b.v. reserves the right to change design and specifications without further notice.

BMS Microscopes b.v.
Essebaan 50
2908 LK Capelle aan den IJssel
The Netherlands
Tel.: +31(0)10 458 42 22
E-mail: sales@bmsmicroscopes.com

Environment protection first!

Your appliance contains valuable materials which can be recovered or recycled. Leave it at a local civic waste collection point.

General troubleshooting tables

If difficulties should be encountered in the course of operation, and no major instrument malfunctions can initially be detected, please recheck the symptoms, referring to the tables provided below, before contacting your service representative.

OPTICAL SYMPTOMS	CAUSES	COUNTERMEASURES
Darkness at the periphery, or uneven viewfield brightness	Revolving nosepiece not in click-stop position (objective not centered in optical path)	Revolve to click-stop position (swing the objective correctly into the optical path)
	Field diaphragm not centered	Centering
	Field diaphragm closed too far	Open properly
	Dirt or dust on the lens(condenser, objective, eyepiece etc)	Cleaning
Dirt or dust in the viewfield	Dirt or dust on the lens (condenser, objective, eyepiece etc)	Cleaning
	Dirt or dust on the slide	Cleaning
	Condenser position too low	Correct positioning
Poor image quality (low resolution or contrast)	No coverglass attached to the slide	Attach coverglass
	Coverglass too thick or thin	Use coverglass of specified thickness (0.17mm)
	Slide upside-down	Turn over the slide
	Immersion oil on dry system objective (especially 40X)	Cleaning
	Air bubbles in immersion oil	Remove bubbles
	Condenser aperture and field diaphragm opened too far	Close properly
	Dirt or dust on the entrance lens	Cleaning
	Condenser aperture closed too far	Open properly
	Condenser position too low	Raise to the position where the diaphragm image is in focus
Image dark on one side	Revolving nosepiece not in click-stop position	Revolve to click-stop position
	Floating specimen	Fasten securely
Image shifts during focus	Specimen rises from stage surface	Place it stable
	Revolving nosepiece not in click-stop position	Revolve to click-stop position
Image tinged yellow	Daylight filter not used	Use daylight filter
Insufficient illumination brightness	Condenser aperture too small	Readjust aperture
	Condenser position too low	Correct positioning
	Dirt or dust on the lens (condenser, objective, eyepiece etc)	Cleaning

MECHANICAL SYMPTOMS	CAUSES	COUNTERMEASURES
Specimen image unsmooth	Mechanical stage not securely fastened	Tighten all fasteners
Image not focusable with high power objectives	Slide upside-down	Turn slide over
	Coverglass too thick	Use coverglass of specified thickness (0.17mm)
Binocular image not integrated	Interpupillary distance not correctly adjusted	Adjustment
Excessive eye fatigue	Incorrect diopter adjustment	Correct adjustment
	Inadequate brightness or illumination	Adjust brightness with control dial

ELECTRICAL SYMPTOMS	CAUSES	COUNTERMEASURES
Lamp does not light when switched ON	No electrical power	Check power cord connection
	LED Lamp not inserted	Insert correctly
	LED failure	Replacement: if the same symptom occurs after replacing the lamp, contact your dealer
Insufficient illumination brightness	Unspecified lamp bulb used	Replace with specified lamp bulb
	Voltage too low	Increase brightness with control dial
Unstable lamp brightness	LED Lamp about to fail	Replacement
	LED Lamp not correctly inserted into socket	Check for positive connection

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