

---

# High Definition Video Microscope

## Manual of Use

### Technical Features

Specification	Model	<b>756460</b> <b>Smartscope</b>
	Series	Measurement HDMI Industrial Camera
	Sensor size	SONY IMX335 sensor, 1/3 Inch
	Effective pixel	Full HD
	Frame rate	1920*1080P@ 60FPS
	Output	HDMI , USB2.0 Type-A, minicard slot
	Pixel size	2.0μm*2.0μm
	Power	DC-12V/2A
Features	Mark function	Point coordinates, Reticle, Coordinates, Text annotation
	Length measurement	Straight length, broken line length, curve line length, distance between parallel lines, point-line distance
	Center distance of circles	Radius setting circle, two points setting circle, three points setting circle
	Geometric measurement	Line length, radius setting circle, two points setting circle, three points setting circle, concentric circles
	Geometry area	Polygon, rectangle
	Image capture	Supports static photos
	Video record	Support 1080P@30 video
	Network video	no
Other features	Horizontal/Vertical mirror, Flip, Monochrome color, Image freeze, Image contrast	
Package	Accessory	Mouse, DC-12V Power adapter, HDMI cable

---

# Production Picture



# Function menu

---


## Boot Screen

Connect the power adapter (12V/2A) , the camera starts up and displays the boot screen as


belows:

## Main features


User move the mouse to left or right after picture display, the main interface pops up(as picture)



## Assistive tools




## Image tools




The icons are in order as above: Zoom in, zoom out, Mirror, Flip, Freeze, Monochrome, Contrast(Click to compare the image from external storage facility with the live image of video)、 Four images contrast: Click the button, and four image will display, double click one of images, the image will be frozen. Double click the image again, unfreeze the image. This helps to compare the products from different angle.

## Parameter setting




Move the mouse to drag the corresponding sidebar under this menu, it can adjust the parameters of the image. If there is a deviation on automatic exposure or automatic white balance, the user can adjust target brightness and exposure time manually to adjust image brightness. At the time of light stability, click one push to adjust the white balance, which can not do auto white balance after calibration no longer

## Image setting



Under this menu, you can drag the progress bar by the mouse to adjust the parameters of relevant images

### System setting



Under this menu, user can set up language, time. View the information of external storage device, camera model, firmware and version. (For time adjustment, the mouse can move to the corresponding position and scroll the mouse wheel to set the corresponding time.) The left and right arrows represent the display position of the switch menu.

## Shortcut


A. Right-click the photo button under the main menu by the mouse, user can set up the camera timer function, and set up the corresponding interval time and quantity of photos.



B. Click the mouse wheel under the main menu to set the photo format (Snap the photo by BMP format, it takes longer time, MJP format is system default)











C. Click the button **mm/px** in the measurement interface, turn on or off the image scale function





D. In the preview state, user scroll the mouse wheel to realize the image digital zoom in and out function quickly.



## Calibration setting

Calibration		
	200mm:200px	   
Item	  	
Name	200mm:200px	
Length	200.000	mm 
Pixel	200.000	px
Scale	1.000000	mm/px

Click the button  to set up the calibration, the user can choose the line and circle calibration alternatively. Input the calibration name in the corresponding name blank, on the preview interface to draw a standard length of the line or circle according to the actual length, and input the number in the length column, then confirm and select the corresponding unit, Click the save button again to complete calibration adding.

Calibration		
	200mm:200px	   



---

## Grid line setting

### Grid

Group  ▼ Line  ▼

Color  ▼ Line  ▼

---

### Measurement attribute

Line c  ▼ Line w 1

Font c  ▼ Font s 16

Click the "All" button to turn on 8 horizontal and vertical scale lines. User can turn on or off any grid line base on your need. User can also set up the color and thickness of each line. Under the "Group" option, it can save 8groups of grid lines, which is convenient to open the set grid lines quickly when testing the product.

In the measurement property, user can select the width and color of the graph line, select the size of the font color of the measurement information.



Microscopes