User manual
Benutzerhandbuch
Manual de usuario
Gebruikershandleiding

Microscopes

Article codes | Artikelnummern | Códigos artículos | Artikelnummers:

74920, 74945, 749455, 74940, 749555 & 749575

WWW.BMSMICROSCOPES.COM


BMS_XL-Stereo_series_GB-DE-ES-NL_manual_v1.PDF9001749

BMS 130, LED
BMS 133 Zoom Trino, LED
BMS 133 Zoom Trino, LED, wireless
BMS 133 Zoom, LED
BMS 140 Zoom, LED
BMS 143 Zoom Trino, LED

GB

Thank you for choosing BMS Microscopes.

The XL stereo microscope, with its versatility in magnification options, adjustable lighting and ergonomic design, is a valuable microscope suitable for various applications. Whether you require intricate detail examination or routine inspections, this microscope proves its worth in fields such as biology, electronics, research and quality control. Its adaptability and high-quality optics make it equally well-suited for scientific research and industrial quality control. Ensuring dependable performance and ergonomic design across diverse fields.

Including: detachable power cord, dust cover, eyecaps, frosted glass disc and manual

Package dimensions (Width x Height x Depth): 29 x 51 x 36 cm

Gross weight: 5,3 kg

Warranty: 5 years

Article code: 74920


STEREO MICROSCOPE BMS 133 ZOOM TRINO, LED, WIRELESS (article code 749455)

Only this special WIRELESS version with article code 749455 can also operate on batteries (both rechargeable and non-rechargeable).

Rechargeable batteries (AA, 1.5V), available as an optional item, can be ordered as a pack of three (item no. 29554). When this battery pack is inserted, it will automatically charge when the microscope is connected to the mains and the power switch of the microscope is turned on. You will know the battery pack is fully charged when the LED indicator on the back of the microscope turns green. The microscope can then be used disconnected from the mains for up to 10 hours. If the LED indicator turns red, it means the batteries need to be recharged. The charging time is approximately 3 hours. The microscope can also operate using non-rechargeable batteries (AA, 1.5V, 3 pieces). PLEASE NOTE: do not connect the microscope to the main electricity when these non-rechargeable batteries are in use, as doing so will cause severe damage to the microscope

BMS 133 Zoom Trino, LED
Article code: 74945

Including: detachable power, cord dust cover, eyecaps, frosted glass and manual

Package dimensions (Width x Height x Depth): 31 x 51 x 37 cm

Gross weight: 5,9 kg (74940) and 6,1 kg (74945 and 749455)

Warranty: 5 years

Article codes: 74940, 74945 and 749455

BMS 143 Zoom Trino, LED
Article code: 749575

Including: dust cover, eyecaps, Ø 95 mm frosted glass, black/white plate and manual

Package dimensions (Width x Height x Depth): 37 x 50 x 27 cm

Gross weight: 6,8 kg (749555) and 7,1 kg (749575)

Warranty: 5 years

Article codes: 749555 and 749575

GETTING STARTED

- **Use the glass stage:** press the glass stage on the sunken place then the other side of the glass stage will be lifted.
- **Set the specimen slide:** set the specimen in the center of stage plate. If necessary, clamp the slide with the clips. Turn on the light.
- **Diopter and focus adjustment:** the focusing knob and observe the specimen through the right eyepiece till the image of the specimen is clear. Observe the specimen through the left eyepiece and adjust the diopter adjustment ring till the image is clear.
- **Adjust the interpupillary distance:** adjust the prism housing along the direction of arrow till the observation is comfortable.
- **Use eyepiece shields:** for user who does not wear glasses, hold the diopter-adjusting ring to prevent them from rotating and turn the eyepiece till the eyepiece shield fit the observer well. For user who wears glasses, take the eyepiece shields off before observation.
- **Installing and removing the optional eyepiece micrometer:** turn and remove the mounting ring from the eyepiece. Clean the eyepiece micrometer and mount it to the mounting ring with the inscription side downward. Gently twist the mounting ring with the eyepiece micrometer into the eyepiece till tightening securely. To remove the eyepiece micrometer, take down the mounting ring by twisting and take out of the micrometer and the wrap it in clean soft paper for storage.
- **Install the illumination device:** insert the illumination device in the bracket with the protrudent side toward the lock-screw and tighten the lock-screw. Put the plug into the socket of the pillar stand.
- **Choose the optical system:** you can alternate the binocular and video capture by pushing or pulling "the pole". You can attain binocular observation by pushing "the pole" inside, or attain video capture by pulling it outside.
- **Adjust the brightness of the bottom light:** turn knob clockwise to increase, counterclockwise to decrease, following base markings.
- **Mount the photo eyepiece and the PK-mount adapter:** put the photo eyepiece into the eyepiece socket of the tri-ocular. Connect the PK-mount adapter with the photo eyepiece and then tighten the lock-screw.
- **Replace the lamps:** press the stage on the sunken place then the other side will be lifted. Take the lamp out of the jack. Put a new lamp into the jack thoroughly. Recover the stage plate. Note: before replacing the lamps, turn off the power first. Avoid violence while the lamp is plugged into the jack.
- **Replace the fuse:** screw the fuse tube out with a screwdriver and pull the fuse out of the tube. Remove the fuse and mount it in an adverse way.

POWER OPTIONS

Main power cable (100~230V/50Hz), available for EU (item no. 76995), UK (item no. 76999), CH (item no. 76998) and USA (item no. 76996).

DOUBLE SAFETY IS GUARANTEED

- Full separation from the main power through DPDT* switch.
- Indicator LED to show if main power is switched on and if optional battery is charged
- Potentiometer for regulation of the High Power LED
- High Power LED for bright image

*DPDT (double pole, double throw) when switched off, the power is completely shut off. The BMS EduLed microscope has a LED illumination light source. This light source is characterized by low maintenance and low power consumption. This LED source has an average life span of 50,000 hours. Besides that, this LED illumination system lacks heat development and is therefore safe to use by children. It is recommended though that an adult is supervising at all time

SAFETY INFORMATION

- Always disconnect mains plug from power point before maintenance

- Work in a clean and proper environment
- Think of the safety of yourself and those working around you
- Clean all parts coming in contact with the observed specimen to prevent infection
- Place the microscope on a sturdy, level table or working surface
- Always use original accessories

CAUTION

The delicate components may get damaged if the microscope is not handled with care. Use the microscope in a clean environment. Do not use it when it is subject to direct sunlight, high temperatures, humidity or dust. Carry the microscope with care, do not hold it by the stage. Remove specimen before transport.

CARE & MAINTENANCE

- Please use our special microscope maintenance kit (item no. 76026).
- Always keep the optical parts clean. Please use the dust cover to protect it.
- Avoid contamination of the objectives.
- In case there is dirt on the lens surface, use distilled water to moist and soft cloth to clean the lens surface. Be very careful and make sure to avoid scratches on a lens surface. Cleaning can also be done by using an alcohol solution or an optical cleaning solution. Ventilate the room. Do not soak the lens in (cleaning) liquid.
- To clean the eyepiece- and objective-lens surface, move the cloth in a spiral motion from the center to the rim.
- Never make a zig-zag movement since this will only spread the dirt.
- To clean painted parts, use also a slightly moistened soft cloth.
- Never use an organic solvent. Use a mild detergent, like soft hand soap.
- Do not remove the objectives from the nosepiece. The objectives were placed in a particular sequence and are adjusted in order to have a parcentered optical system.
- Never dismantle the microscope. This can only be done by factory and/or dealer authorized staff.

SPECIFICATIONS

For more specifications, please view the product page on the website: www.bmsmicroscopes.com

It is within our company policy to develop continuously. BMS Microscopes b.v. reserves the right to change design and specifications without further notice.

BMS Microscopes b.v.

Essebaan 50

2908 LK Capelle aan den IJssel


The Netherlands

Tel.: +31(0)10 458 42 22

E-mail: sales@bmsmicroscopes.com

Environment protection first!

Your appliance contains valuable materials which can be recovered or recycled. Leave it at a local civic waste collection point.

General troubleshooting tables

If difficulties should be encountered in the course of operation, and no major instrument malfunctions can initially be detected, please recheck the symptoms, referring to the tables provided below, before contacting your service representative.

SYMPTOMS	CAUSES	COUNTERMEASURES
Double images	Interpupillary distance is not correct	Readjust it
	Diopter adjustment is not correct	Readjust it
	Magnification of each eyepiece is not the same size	Mount the same size eyepiece
Dirt appears in the field of view	Dirt on the specimen	Clean the specimen
	Dirt on the surface of eyepiece	Clean the surface
Image is not clear	Dirt on the surface of the objective	Clean the objective
Image is not clear while the focus changing	Diopter adjustment is not correct	Readjust the diopter
	Focus is not correct	Readjust the focus
The focusing knob is not smooth	The focusing knob is too tight	Loosen it to a suitable position
The image is obscure because of the head slipping down by itself during observation	The focusing knob is too loose	Tighten it to a suitable position
Incision image appears in the field of view or of the video view	The pole is not in correct position	Pull or push it to the correct position
Eyes fell tired easily	Diopter adjustment is not correct	Adjust the diopter
	Brightness of light is not correct	Adjust the brightness
Bulb does not work when the switch is on	No power on	Check the connection with the power supply
	The bulb was not inserted correct	Insert it correctly
	Bulb is wrong	Replace with a new one
Bulb is burned out suddenly	Use the wrong bulb	Replace with a correct one
	The voltage is too high	Control the voltage: e.g. use voltage regulator
Brightness is not enough	Use the wrong bulb	Replace with a correct one
	The voltage is too low	Increase the input voltage
The bulb flickers or the brightness is unstable	The bulb will burn out soon	Replace with a new one
	The bulb was not inserted correctly	Insert it correctly

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