User manual
Benutzerhandbuch
Manual de usuario
Gebruikershandleiding

Microscopes

Article codes | Artikelnummern | Códigos artículos | Artikelnummers:
8746110, 8746110-60, 8746110-100, 8746120, 8746120-60, 8746120-100,
8746130, 8746130-60, 8746130-100, 9746120, 9746140 & 9746130

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
BMS_D2-D3_series_GB-DE-ES-NL_manual_v1.PDF9001D2D3

GB

Thank you for choosing BMS Microscopes.

The BMS D2 and D3 microscope make it possible to study structures and textures of botanical and zoological nature at cell level. Both prepared slides or self-made slides can be observed. The material used should allow light to pass easily, using a specimen that is too thick will result in a blurred and undefined image. This microscope is developed according to the latest technology on optical, electrical and mechanical details. The BMS D2 and D3 are equipped with an inventive power system.

Including: dust cover and manual

Package dimensions (Width x Height x Depth) 41 x 52 x 27 cm

Gross weight: 7,2 kg (mono) and 7,8 kg (bino/trino)

Warranty: 5 years

Article codes: 8746110, 8746110-60, 8746110-100, 8746120, 8746120-60, 8746120-100, 8746130, 8746130-60 and 8746130-100


BMS D2
MONO

BMS D2 or D3
BINO

BMS D2 or D3
TRINO

Eyepiece: WFP 10x/22 mm with foldable eyecups

Tube: 45° inclined

Head: bino or trino, 360° rotatable

Nosepiece: quintuple, ball bearing type with rubber grip

Stage: Rackless, 144 x 200 mm, suitable for two slides at one time.

Equipped with integrated movable specimen holder with ergonomically placed coaxial controls for accurate X-Y movement.

Focusing: coaxial coarse and fine adjustment on both sides

Anti-fungus optics

Objectives: DIN ePlan High Resolution Infinity. Magnification: 40x, 100x and 400x. Optional to be extended to 2000x.

Condenser: Abbe condenser (N.A. 1.25) in holder, with centering screws. Adjustable height via rack and pinion.

New generation High Power LED. Lifespan LED: approx. 50.000 hours.

Anti-bacterial paint

Including: dust cover and manual

Package dimensions (Width x Height x Depth): 46 x 56 x 26 cm

Gross weight: 8,7 kg

Warranty: 5 years

Article codes: 9746120, 9746140 and 9746130

GETTING STARTED

- LED-lamp illumination: turn the brightness control dial forward to achieve the desired brightness.
- Specimen mounting: mounting with the cover glass facing upwards.
- Focusing with 10X objectives: choose the objective 4X or 10X. Turn the coarse focusing knob to get the objective focused on the slide and then turn the fine focus adjustment knob to get the image sharp and clear.
- Interpupillary distance adjustment: adjust the interpupillary distance of binocular until the right-left field of view can be composed one.
- Eyepiece diopter adjustment: firstly, observe at right tube with right eye, adjust coarse-fine focusing knob to image clearly. Secondly, observe at left tube with left eye, adjust the diopter control 1 to image clearly.
- Objective selection: rotate the revolving nosepiece to the desired objective magnification, assuring it click-stops into position.
- Aperture diaphragm adjustment: this diaphragm is designed for numerical aperture adjustment, not brightness. Generally, a good image of sufficient contrast can be achieved with the aperture diaphragm stopped down to approximately 70-80% of the objective's numerical aperture. To observe the diaphragm image, remove the eyepiece and look down the open tube at the exit pupil of the objective.
- Field diaphragm adjustment: for the field lens unit with field diaphragm; stop down the field diaphragm until the illuminated area closely circumscribes the eyepiece viewfield. The field diaphragm is used to control the specimen's illuminated area relative to the microscope viewfield. If it is opened to a larger aperture than necessary, extraneous light will enter the viewfield, reducing image quality and contrast.

POWER OPTIONS

1. Main power cable (100~230V/50Hz), available for EU (item no. 76995), UK (item no. 76999), CH (item no. 76998) and USA (item no. 76996).
2. Rechargeable batteries (AA, 1.5V), optional item, can be ordered as a pack of three (item no. 29554). When this battery pack is placed, it will be charged automatically when the microscope is connected to the mains and the power switch of the microscope is switched on. When the LED indicator, located on the back of the microscope, is green, the battery pack is fully charged. The microscope can now be used, disconnected from the mains, for 10 hours. When the LED indicator is red, the batteries need to be charged again. Charging time is about 3 hours.
3. Non rechargeable batteries (AA 1.5V, 3 pieces). The microscope will also work with these batteries. WARNING: do not connect the microscope to the main electricity when these batteries are placed. This will cause severe damage to the microscope.

DOUBLE SAFETY IS GUARANTEED

- Full separation from the main power through DPDT* switch.
 - Indicator LED to show if main power is switched on and if optional battery is charged
 - Potentiometer for regulation of the High Power LED
 - High Power LED for bright image
- *DPDT (double pole, double throw) when switched off, the power is completely shut off. The BMS D2 and D3 microscopes have a LED illumination light source. This light source is characterized by low maintenance and low power consumption. This LED source has an average life span of 50,000 hours. Besides that, this LED illumination system lacks heat development and is therefore safe to use by children. It is recommended though that an adult is supervising at all time.

SAFETY INFORMATION

- Always disconnect mains plug from power point before maintenance
- Work in a clean and proper environment
- Think of the safety of yourself and those working around you
- Clean all parts coming in contact with the observed specimen to prevent infection

- Place the microscope on a sturdy, level table or working surface
- Always use original accessories

CAUTION

The delicate components may get damaged if the microscope is not handled with care. Use the microscope in a clean environment. Do not use it when it is subject to direct sunlight, high temperatures, humidity or dust. Carry the microscope with care, do not hold it by the stage. Remove specimen before transport.

REPLACING FUSE

- Always disconnect mains plug from power point before any kind of maintenance. The BMS D2 and D3 microscopes are equipped with one glass fuse, 1,0 A.
- Wait until the microscope and electrical parts have cooled down.
 - The fuse is located on the bottom of the microscope
 - Open the housing of the fuse by turning it counter clockwise.
 - Take broken fuse out.
 - Replace new fuse after being sure of using same rating in new fuse.
 - Screw fuse holder clockwise back in place.

CARE & MAINTENANCE

- Please use our special microscope maintenance kit (item no. 76026).
- Always keep the optical parts clean. Please use the dust cover to protect it.
- Avoid contamination of the objectives.
- In case there is dirt on the lens surface, use distilled water to moist and soft cloth to clean the lens surface. Be very careful and make sure to avoid scratches on a lens surface. Cleaning can also be done by using an alcohol solution or an optical cleaning solution. Ventilate the room. Do not soak the lens in (cleaning) liquid.
- To clean the eyepiece- and objective-lens surface, move the cloth in a spiral motion from the center to the rim.
- Never make a zig-zag movement since this will only spread the dirt.
- To clean painted parts, use also a slightly moistened soft cloth.
- Never use an organic solvent. Use a mild detergent, like soft hand soap.
- Do not remove the objectives from the nosepiece. The objectives were placed in a particular sequence and are adjusted in order to have a parcentered optical system.
- Never dismantle the microscope. This can only be done by factory and/or dealer authorized staff.

SPECIFICATIONS


For more specifications, please view the product page on the website: www.bmsmicroscopes.com

It is within our company policy to develop continuously. BMS Microscopes b.v. reserves the right to change design and specifications without further notice.

BMS Microscopes b.v.
Essebaan 50
2908 LK Capelle aan den IJssel
The Netherlands
Tel.: +31(0)10 458 42 22
E-mail: sales@bmsmicroscopes.com

Environment protection first!

Your appliance contains valuable materials which can be recovered or recycled. Leave it at a local civic waste collection point.

General troubleshooting tables

If difficulties should be encountered in the course of operation, and no major instrument malfunctions can initially be detected, please recheck the symptoms, referring to the tables provided below, before contacting your service representative.

OPTICAL SYMPTOMS	CAUSES	COUNTERMEASURES
Darkness at the periphery, or uneven viewfield brightness	Revolving nosepiece not in click-stop position (objective not centered in optical path)	Revolve to click-stop position (swing the objective correctly into the optical path)
	Field diaphragm not centered	Centering
	Field diaphragm closed too far	Open properly
Dirt or dust in the viewfield	Dirt or dust on the lens (condenser, objective, eyepiece etc)	Cleaning
	Dirt or dust on the slide	Cleaning
	Condenser position too low	Correct positioning
Poor image quality (low resolution or contrast)	No coverglass attached to the slide	Attach coverglass
	Coverglass too thick or thin	Use coverglass of specified thickness (0.17mm)
	Slide upside-down	Turn over the slide
	Immersion oil on dry system objective (especially 40X)	Cleaning
	Air bubbles in immersion oil	Remove bubbles
	Condenser aperture and field diaphragm opened too far	Close properly
	Dirt or dust on the entrance lens	Cleaning
	Condenser aperture closed too far	Open properly
Image dark on one side	Condenser position too low	Raise to the position where the diaphragm image is in focus
	Revolving nosepiece not in click-stop position	Revolve to click-stop position
	Floating specimen	Fasten securely
Image shifts during focus	Specimen rises from stage surface	Place it stable
	Revolving nosepiece not in click-stop position	Revolve to click-stop position
Image tinged yellow	Daylight filter not used	Use daylight filter
Insufficient illumination brightness	Condenser aperture too small	Readjust aperture
	Condenser position too low	Correct positioning
	Dirt or dust on the lens (condenser, objective, eyepiece etc)	Cleaning

MECHANICAL SYMPTOMS	CAUSES	COUNTERMEASURES
Specimen image unsmooth	Mechanical stage not securely fastened	Tighten all fasteners
Image not focusable with high power objectives	Slide upside-down	Turn slide over
	Coverglass too thick	Use coverglass of specified thickness (0.17mm)
Binocular image not integrated	Interpupillary distance not correctly adjusted	Adjustment
Excessive eye fatigue	Incorrect diopter adjustment	Correct adjustment
	Inadequate brightness or illumination	Adjust brightness with control dial

ELECTRICAL SYMPTOMS	CAUSES	COUNTERMEASURES
Lamp does not light when switched ON	No electrical power	Check power cord connection
	LED Lamp not inserted	Insert correctly
	LED failure	Replacement: if the same symptom occurs after replacing the lamp, contact your dealer
Insufficient illumination brightness	Unspecified lamp bulb used	Replace with specified lamp bulb
	Voltage too low	Increase brightness with control dial
Unstable lamp brightness	LED Lamp about to fail	Replacement
	LED Lamp not correctly inserted into socket	Check for positive connection

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