

# BMS DISSECTION



**BMS**

Microscopes

## PURE WHITE BMS COVER SLIPS


### Article code 11290

- 0.13 - 0.17 mm
- 18 x 18 mm
- 100 pcs per box



### Article code 11291

- 0.13 - 0.17 mm
- 20 x 20 mm
- 100 pcs per box




### Article code 11292

- 0.13 - 0.17 mm
- 22 x 22 mm
- 100 pcs per box



### Article code 11397

- 0.13 - 0.17 mm
- Ø 16 mm
- 100 pcs per box



### Article code 12288

- Thickness: 1.0 - 1.2 mm
- Size: 26 x 76 mm
- Ground edges
- With fine frosted marking area on both sides
- ISO 9001
- 50 pcs per box



### Article code 12281

- Thickness: 1.0 - 1.2 mm
- Size: 76 x 26 mm
- Ground edges
- Plain without frosted end
- ISO 9001
- 50 pcs per box



# BMS SLIDES & COVER SLIPS



## Article code 12283

- Thickness: 1.0 - 1.2 mm
- Size: 76 x 26 mm
- Ground edges
- With fine frosted marking area
- ISO 9001
- 50 pcs per box

## Article code 12290

- Thickness: 1.0 - 1.2 mm
- Size: 76 x 26 mm
- Ground edges
- With ground-in and polished cavity Ø 15 - 18 mm
- ISO 9001
- 50 pcs per box



## Article code 12286

- Thickness: 1.2 - 1.4 mm
- Size: 75 x 26 mm
- Ground edges
- Plain without frosted end
- ISO 9001
- 50 pcs per box



## Article code 12287

- Thickness: 1.2 - 1.4 mm
- Size: 76 x 26 mm
- Ground edges
- With fine frosted marking area
- ISO 9001
- 50 pcs per box

## Article code 12289

- Thickness: 1.2 - 1.4 mm
- Size: 75 x 26 mm
- Ground edges
- With fine frosted marking area on both sides
- ISO 9001
- 50 pcs per box



## BMS DISSECTING CASES



### Dissecting set in case, 5-parts (article code 166455)

- Scissors, straight, sharp, 105 mm
- Forceps, straight, fine tip, 130 mm
- Scalpel handle no.4, 110 mm
- Needle with holder, straight
- Scalpel no. 23, lancet-shaped, 5 pieces



### Dissecting set in case, 9-parts (article code 16642)

- Needle, straight (2x)
- Needle, lancet-shaped
- Forceps, straight, blunt, 130 mm
- Forceps, straight, fine tip, 130 mm
- Scissors, straight, sharp points, 105 mm
- Scalpel handle no.4, 110 mm
- Scalpel No. 23, lanced-shaped, 5 pieces
- Spatula
- Dissecting case



### Dissecting set in case, 9-parts (article code 16645)

- Forceps, straight, fine tip, 130 mm
- Scissors, straight, sharp points, 105 mm
- Scissors, curved, sharp points, 105 mm
- Sectioning razor
- Scalpel handle no.4, 110 mm
- Needle, straight (2x)
- Needle, lancet-shaped
- Scalpel No. 23, lancet-shaped, 5 pieces
- Dissecting case



### Dissecting case (article code 16644)

- Made of synthetic leather)
- 20 x 31 cm (open)



# BMS MAGNIFIERS



Aplanatic hand magnifier  
(article code 15907)

- Lens Ø: 18 mm
- Magnification: 10x



Folding pocket magnifier one  
lens (article code 15890)

- Lens: Ø 22 mm
- Magnification: 4x



Folding pocket magnifier two  
lenses (article code 15891)

- Lens: Ø 22 mm
- Magnification: 4x & 8x



Folding pocket magnifier three  
lenses (article code 158925)

- Lens: Ø 22 mm
- Magnification: 5x, 10x & 15x



Folding pocket magnifier triplet  
(article code 15896)

- Lens: Ø 18 mm
- Magnification: 10x




Hand magnifier 3x/ Ø 50 mm  
(article code 15904)

- Lens Ø: 50 mm
- Magnification: 3x



Hand magnifier 3x/ Ø 75 mm  
(article code 15905)

- Lens Ø: 75 mm
- Magnification: 3x



Folding linen tester magnifier  
(article code 15912)

- Lens Ø: 19 mm
- Magnification: 7.5x



Stand magnifier 2.5x and 6x/ Ø  
90 mm (article code 15900)

- 3-way head
- Lens Ø: 90 mm
- Magnification: 2.5x and 6x
- With LED illumination and  
batteries

# BMS FORCEPS & NEEDLES

	STRAIGHT & BLUNT	STRAIGHT & FINE TIP	CURVED & FINE TIP
105 mm	16580	16590	
130 mm	16581	16591	16596
145 mm		16592	16597
160 mm	16582		

	50 mm	35 mm	40 mm
	NO HANDLE	PLASTIC HANDLE	STAINLESS STEEL HANDLE
STRAIGHT	16625	16610	16609
LANCET-SHAPED	16626	16615	16614



Article code 16598

Forceps, surgical/  
tissue, 140 mm



Article code 16605

Forceps, splinter,  
140 mm



Article code 16600

Forceps, cover-  
glass, 105 mm



Article code 16606

Forceps, galvanic,  
90 mm



Article code 16627

Inoculating loop

- Length: 85 mm
- Ø: 0,5 mm



Article code 16620

Needle holder

For 16625, 16626  
and 16627



Article code 16611

Seeker/Probe with  
handle

- Length: 145



Article code 16630

Double End  
seeker/probe

- Length: 130



Article code 16884

- Straight
- Sharp points
- 105 mm




Article code 16885

- Curved
- Sharp points
- 105 mm



Article code 16883

- Sharp
- Blunt points
- 140 mm



Article code 16886

- Scissors paper
- 170 mm

# BMS DISSECTION



Article code 16629  
Cartilage knife



Article code 16905  
Scalpel stainless



Article code 16911  
Scalpel no. 23

- Steel
- Lancet-shaped
- 5 pieces
- For 16910 no.4



Article code 16912  
Scalpel no. 21/24

- Steel
- Curved
- 5 pieces
- For 16910 no. 4



Article code 16917  
Razor blades with security strip

- 5 pieces



Article code 16646  
Dissecting tray with wax

- Work surface: 13,5 x 23,7 cm
- Type: beeswax
- Material tray: stainless steel
- Product dimensions: 265 x 162 x 20 mm
- Well suited for laboratory work and school purposes
- Reusable tray after melting the wax




Article code 16909  
Scalpel handle no. 3

- Size: 100 mm



Article code 16910  
Scalpel handle no. 4

- Size: 110 mm



Article code 16616  
Insect pins, 100 pieces, no. 00

- Stainless steel
- Length: 38 mm
- Thickness: 0,30 mm



Article code 16618  
Insect pins, 100 pieces, no. 4

- Stainless steel
- Length: 38 mm
- Thickness: 0,55 mm



Article code 16635  
Sectioning Razor

Flat one side, the other being hollow ground



Article code 16648  
Handmicrotome, cylindrical

Handmicrotome with table clamp and blade



# BMS BUG VIEWERS




Magnifying bug viewer  
(article code 778165)

- Magnification: 3x
- Size Ø: 71
- Grid: yes
- Two sections
- Plastic spider



Bug viewer 2x/4x, 150 ml  
(article code 778125)

- Magnification: 2x and 4x
- Width by Height: 71 x 75 mm
- Bi-convex: 47 mm
- Cup height: 67 mm
- Three sections
- Plastic spider



Magnifier 3 in 1  
(article code 74212)

- Sizes Ø: 100, 20 and 20 mm
- Magnification: 2x, 3x & 4x
- Ø 100 mm: magnification 2x
- Ø 20 mm: magnification 3x
- Ø 20 mm: magnification 4x



Magnifier 3x and 6x  
(article code 15903)

- 2 in 1 magnifier
- Sizes Ø: 60 and 20 mm
- Magnification: 3x and 6x



Three section bug viewer, retractable  
(article code 742155)

- Magnification: 2,5x and 4,5x
- Lens: fixed lens Ø 47 mm and flip lens of Ø 41 mm
- Sizes: 70x35 mm (WxH)
- Plastic spider



Magnifying glass on stand  
(article code 15902)

- Hands-free observation
- Magnification: 5x
- Housing Ø: 100 mm
- Size Ø: 75 mm



Two way bug viewer  
(article code 77807)

- Magnification 4x
- Top and bottom view
- Easy grip handles
- Plastic spider



Microscopes

Essebaan 50  
2908 LK Capelle aan den IJssel  
THE NETHERLANDS

+31 (0)10 45 842 22 | [info@bmsmicroscopes.com](mailto:info@bmsmicroscopes.com) | [www.bmsmicroscopes.com](http://www.bmsmicroscopes.com)